

DermAssist®

Latex Surgical Gloves Series 133

These low-protein, powder-free gloves eliminate airborne particles of latex protein associated with powdered surgical gloves. Ultra-low modulus is designed for greater comfort, less hand fatigue during long surgical procedures.

Characteristics: Sterile, Powder-Free



SIZE	ITEM NO.
5.5	133550
6	133600
6.5	133650
7	133700
7.5	133750
8	133800
8.5	133850
9	133900

All Sizes:
50 Pairs/Box, 4 Boxes/Case



Features:

- Textured Finish = Improved Wet/Dry Grip
- Low Protein = Reduced Changes of Latex Allergic Reactions
- Low Modulus = A Softer, More Comfortable Fit

Manufactured in compliance with multiple international standards, including ASTM D3577, ASTM D5151, ASTM F1671, and ASTM D5712

Tensile Strength & Elongation	Before Aging	After Accelerated Aging
Tensile Strength (Mpa)	30	24
ASTM Requirement Min. (Mpa)	24	18
Elongation (%)	1290	1150
ASTM Requirement Min. (%)	750	560

Average Length	11.5 in ♦ 290 mm
Average Palm Thickness	7.5 mil ♦ 0.19 mm
Average Finger Thickness	7.9 mil ♦ 0.20 mm



Intertek
Innovative Healthcare Corporation is certified to ISO 13485:2016 QMS for medical devices.



Product Solutions You Trust

Canadian Medical Supplies
Tel: 647-514-6290

View this product on our website:

